

Blood Vessel Volumetric Data Pipeline

笔记本: Cornell Univ.

创建时间: 2025/6/30 10:47

更新时间: 2025/10/20 17:17

作者: 刘驰

1. Skeleton Extraction

- Perform 3D median filter on 3D structure stack (radius: 2) ;
- Enter \Step1 - Vessel_skeletonization, Run:

```
GU_extractVesselness_Skel('sampledata\mouse1_fineStack_Hz_740-675um_6pt5mW_00001_med2.tif');
```

- Use 'mouse1_fineStack_Hz_740-675um_6pt5mW_00001_med2_skelInterp.tif' as skeleton.

2. Channel Split

- Enter Vitural Env: conda activate dataprocess

```
python bidirectional_tiff_deinterleaver.py 'E:\AES Data\20250604 practice wild type mouse blood vessel volume imaging\mouse imaging\mouse1_AES_dualPort_50pt4Hz_684-732um_3pt4mW_0pt53mW_00001.tif' 'E:\AES Data\20250604 practice wild type mouse blood vessel volume imaging\mouse imaging\file1' --num-chunks 5
```

- Perform median filtering(radius:1), then Generate MIPs of 6 planes

```
python process_volume_median_mip.py
```

3. Find segment in Skeleton

- Use 'mouse1_fineStack_Hz_740-675um_6pt5mW_00001_med2_skelInterp.tif' as skeleton. Perform temporal-color encode, then record start, end, and possible waypoint of each segment.
- Save Start and end point in Segments.txt
- Find segments with SegmentsExtraction.m
- Segments Prompt saved to SegmentsPrompt.txt

4. Kymograph extraction with \Step4 -

Kymograph_Extraction\kymograph_from_file.py

```
python kymograph_from_file.py
```

5. Nonlinear Transform Extracted Kymograph

- Transform Extracted kymograph with \Step5 - Nonlinear_Transform\nonlinear_transform_of_3d_graph_segment.m
- \Kymograph_Transformed_cut to modify kymo length

6. Calculate blood flow speed

- For multi-segments calculation, LSPIV_parallel_simple_script_batch.m
- Find proper parameters for each segments, LSPIV_parallel_simple_script_singlesegment.m
- Plot speed, LSPIV_parallel_simple_script_speed_plot.m